



LDLS™ Laser-Driven Light Source

From The Innovators in Light

The LDLS Product Range



•EQ-99FC LDLS™ System



•Compact LDLS with fiber-coupled output



•EQ-99 LDLS™ System



•Compact LDLS for free-space optics

•EQ-1500 LDLS™ System




•High brightness LDLS for free-space optics

Today's Lamp Technologies



• Mercury Discharge



• Excimer



• Deuterium



• Halogen



• Xe & Hg Arc



• Metal Halide

All have limitations in:

- Brightness
- Lifetime/Stability
- Spectral Range


Limitations with Arc-Lamps

- Arc-Lamps use electrodes to conduct electrical current through the gas heat the gas to high temperature
 - The electrodes limit the temperature of the gas
 - ❖ Limited temperature leads to limited DUV
 - The electrodes erode and shorten lamp lifetime
 - Arc flicker reduces effective brightness and adds noise





Broadband	✓
Brightness	✓
Stability	✗
Lifetime	✗

• DUV output falls below 300nm



Limitations with Deuterium Lamps





Source: Hamamatsu

Broadband	X
Brightness	X
Stability	✓
Lifetime	X

- Limited wavelength range ~190nm – 400nm
- Low brightness, relatively large plasma
- Short life: 500-1000hrs to 50% output

Limitations with Tungsten-Halogen Lamps



Source: Ocean Optics

Broadband	✓
Brightness	✗
Stability	✓
Lifetime	✗

- Short lifetime, large filament area, low power <400 nm
- Usually combined with Deuterium lamp to cover broad spectrum
 - Gap in spectrum around 400nm
 - Low brightness from 2 separate emitters
 - Changing spectrum from 2 different lifetimes

A Timely Convergence....

- Traditional Lamp sources have reached the end of their roadmap
 - Incremental improvement in recent years

+

- Costs for infra-red laser diodes have fallen dramatically in recent years
 - Driven by telecom and fiber lasers for industrial applications
 - Costs continue to decline substantially each year.

=

- Catalyst for a new class of light sources that disrupts the status quo

The brightest,
most consistent
broadband light sources
since the sun

EQ-99




EQ-99FC



EQ-1500




Laser-Driven Light Source: Principle of Operation



Light from Arc-Lamp and LDLS

•Xe Arc-Lamp



•Large plasma
limits
brightness

•Arc position
Instability

•Note: anode
glowing red

•LDLS




•Small plasma
...high
brightness

•Highly stable
position

•Note: anode
NOT glowing

- High brightness: ~100 um diameter Xenon plasma,
- Efficient coupling into small fibers or spectrometer slits
- Point source enables collimation over long distances

EQ-99 Spectral Distribution




•Typical Data

EQ-99FC Spectral Distribution

EQ-99FC Typical Performance:


with 230 μ m diameter, 0.22NA, 1m long, DUV/Vis solarization resistant fiber



Power ratios	
Fiber (μ m)	Multiplier
230	1
450	2.02
600	2.26


•Typical Data

EQ-1500 Spectral Distribution



•Typical Data

Comparing LDLS with Traditional Lamps




•Spectral radiance calibrated at 254nm, Typical Data

LDLS: UV-Vis-NIR Radiance

170nm – 1700nm


ENERGETIQ



•Typical Data for EQ-99



LDLS™: Stable & Long-Life

ENERGETIQ




<u>Light Source</u>	<u>Change in Broadband Output /1000 Hrs</u> <u>(Typical)</u>	<u>Life-Test Hours</u> <u>to Date</u>	<u>Notes</u>
EQ-1500	~ -1%	>10,000	Test on-going
EQ-99	~ -1%	>6,000	Test on-going
30W D2 Lamp	-50% (depending on model)		Source: Heraeus Data Sheet
75W Xe Lamp	-25% to -50% (depending on model)		Source:Hamamatsu Data Sheet


EQ-99 Plasma Images

	PS37LH34	PS38LH39
FWHM horizontal, μm	64	61
FWHM vertical, μm	147	140
Images		

Spatial Stability of Arc Lamp vs LDLS




450W Xenon Lamp




EQ-1500

Spatial Stability Results



- Collected and stored 2500 images @ 200 frames per second
- Calculated center of mass for each image using ImageJ (image analysis software)
- Standard deviation of the plasma light intensity center of mass position
 - Horizontal: $0.145 \mu\text{m}$
 - Vertical: $0.094 \mu\text{m}$

Measurement of Temporal Variation in Radiant Flux for EQ-99



- Light collected by 0.3NA Achromat pair from 100μm plasma region
- Calibrated integrating sphere/fiber spectrometer:
 - 1000 spectral samples (400-830nm), 8ms integration time, 8 seconds total
- Measurement repeated for ten EQ-99 units.

Temporal Variation in Radiant Flux From 100 μ m Plasma, 0.3NA

200 μ m diameter pinhole with 2X optics,
400nm to 830nm wavelength band,
1000 samples, 8ms integration time, 8s total

EQ-99	Intensity Variation, 3x std.dev. [%]	P-P/Mean (%)
Unit-A	0.27	0.58
Unit-B	0.33	0.73
Unit-C	0.17	0.34
Unit-D	0.20	0.39
Unit-E	0.19	0.42
Unit-F	0.29	0.65
Unit-G	0.34	0.72
Unit-H	0.21	0.48
Unit-I	0.18	0.42
Unit-J	0.33	0.79
Average	0.25%	0.55%



- For the 8 second sample period, average intensity variation of:
 - <0.1% (1σ)
 - 0.25% (3σ)
 - 0.55% (p-p/mean)

Benefits of LDLS Technology

- Very high brightness across complete spectrum
 - 170nm through visible and out to 2100 nm
- Eliminates need for multiple lamps (replaces D2/Tungsten/Xenon Arc)
 - Simplified optical system
- Excellent Spatial stability
 - Repeatable measurements
- Superior short and long term power stability
 - Repeatable measurements
- Electrodeless operation for long life
 - Reduced consumable costs
 - Minimal recalibration of instrument

Broadband	✓
Brightness	✓
Stability	✓
Lifetime	✓

Applications

- UV-Vis Spectrometry
- Monochromator Source
- PEEM
- Atomic Absorption Spectroscopy
- Materials Characterization
- Environmental Analysis
- Hyperspectral Imaging
- Gas Phase Measurements
- Advanced Microscopes
- Endoscopes/Borescopes



•EQ-99 LDLS System



•EQ-1500 LDLS System


NEW! EQ-99 Manager

- A Smart Controller for EQ-99 Series LDLS™ Products



- Easy on/off capabilities
- Lamp hour meter
- Available Shutter
- Computer control

LDLS Fiber Optic Assemblies



- Turn-Key system with EQ-99FC
- For high brightness applications
- Proprietary termination process
- Longer life than standard fibers
- Choice of Fibers
 - DUV/Vis Solarization Resistant
 - ❖ (180nm – 900 nm)
 - ❖ 115 μ m, 230 μ m, 455 μ m Core
 - Broadband Version
 - ❖ (200nm – 2100 nm)
 - ❖ 100 μ m, 200 μ m, 400 μ m, 600 μ m Core


EQ-99CAL Calibrated LDLS



- Based on Award Winning EQ-99 Design
- Broadband, High Brightness Source
- Calibration Wavelengths: 200nm - 800nm (with single source)
- Irradiance Calibration Traceable to NPL
- Temperature Controlled Lamp House for highly stable measurements
- Long Intervals between Recalibration
 - Recommended 1 year or 1000 hours of operation (which ever comes first)

New EQ-99CAL Calibrated LDLS

Typical Calibrated Irradiance Spectrum, 200nm- 800nm




Coupling Accessories



EQ-99 LDLS with OAP System

- Designed to couple EQ-99 and EQ-1500 light into spectrometers or optical fibers
- Off-Axis Parabolic (OAP) Mirror Assembly
- Efficiently couple full wavelength range
- Variety of focal length and NA options
- Free beam or SMA fiber coupling options
- Housings are Nitrogen purgeable



EQ-1500 LDLS with OAP System

Broadband Fiber Collimators For LDLS



EQ-99FC LDLS with Fiber Collimator

- Designed to collimate EQ-99FC fiber output
- Reflective Mirror Assembly
- Efficiently collimate from 200 – 2400 nm
- Two output beam sizes available; 6.7 mm and 11 mm diameter
- SMA fiber connection
- Low divergence (depends on fiber size)



EQ-99-COL-11
EQ-99-COL-6

Summary

- Very high brightness across complete spectrum
 - 170nm through visible and out to 2100 nm
 - Easy coupling to small fibers and spectrometer slits
 - Ease of collimation
- Eliminates need for multiple lamps (replaces D2/Tungsten/Xenon Arc)
 - Simplified optical system
- Excellent Spatial stability
 - Repeatable measurements
- Superior short and long term power stability
 - Repeatable measurements
- Electrodeless operation for long life
 - Reduced consumable costs
 - Minimal recalibration of instrument



LDLS™ Laser-Driven Light Source
www.energetiq.com